Citoqueratina 8/18 – Anticorpo Monoclonal anti-humano

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1. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 73:626 (1980).
2. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 14:767 (1983).
3. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15.eds.
4. Angus B, Purvis J, Stock D et al. NCL-5D3: a new monoclonal antibody recognizing low molecular weight cytokeratins effective for immunohistochemistry using fixed paraffinembedded tissue. *Journal of Pathology*. 153(4):377- 384 (1987).
5. Angus B, Kiberu S, Purvis J et al. Cytokeratins in cervical dysplasia and neoplasia: a comparative study of immunohistochemical staining using monoclonal antibodies NCL-5D3, CAM 5.2, and PKK1. *Journal of Pathology*. 155(1):71-75 (1988).
6. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 7(9). Order code M29-P (1991).
7. Martin de las Mulas J, Espinosa de los Monteros A, Bautista MJ et al. Immunohistochemical distribution pattern of intermediate filament proteins and muscle actin in feline and human mammary carcinomas. *J Comp Pathol*. 111(4):365-381 (1994).
8. Torii H, Nakagawa H and Ishibashi Y. So-called mixed tumor of the skin on the wrist: an immunohistochemical study. *Journal of Dermatology*. 22(3):213-217 (1995).
9. Maurizii MG, Saverino O and Taddei C. Cytokeratin cytoskeleton in the differentiating ovarian follicle of the lizard *Podaris sicula* Raf. *Molecular Reproduction and Development*. 48(4):536-542 (1997).
10. Hinze P, Feyler S, Berndt J et al. Malignant myoepithelioma of the vulva resembling a rhabdoid tumour. *Histopathology*. 35(1):50-54 (1999).
11. Mehes G, Witt A, Kubista E et al. Circulating breast cancer cells are frequently apoptotic. *American Journal of Pathology*. 159(1):17-20 (2001).
12. Martin CA, Salomoni PD and Badran AF. Cytokeratin immunoreactivity in mouse tissues: study of different antibodies with a new detection system. *Appl Immunohistochem Mol Morphol*. 9(1):70-73(2001).