## SPARC (Osteonectina) – Anticorpo Policlonal anti-humano

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1. Swaroop A, Hogan BL, Francke U. Molecular analysis of the cDNA for human SPARC/osteonectin/BM-40: sequence, expression, and localization of the gene to chromosome 5q31-q33. *Genomics*. 1988;2:37-47
2. Schulz A, Jundt G, Berghäuser KH, Gehron-Robey P, Termine JD. Immunohistochemical study of osteonectin in various types of osteosarcoma. *American Journal of Pathology*. 132(2): 233-238. 1988.
3. Mundlos S, Schwahn B, Reichert T, Zabel B. Distribution of osteonectin mRNA and protein during human embryonic and fetal development. *The Journal of Histochemistry and Cytochemistry*. 40(2): 283-291. 1992.
4. Frizell E, Liu SL, Abraham A, Ozaki I, Eghbali M, Sage EH, Zern MA. Expression of SPARC in normal and fibrotic livers. *Hepatology*. 21 (3): 847-854. 1995.
5. Peggy L. Porter PI, Sage EH, Lane TE, Funk SE, Gown AM. Distribution of SPARC in Normal and Neoplastic Human Tissue. *The Journal of Histochemistry and Cytochemistry*. 43(8): 791-800. 1995.
6. Fanburg JC, Rosenberg AE, Weaver DL, Leslie KO, Mann KG, Taatjes DJ, Tracy RP. Osteocalcin and osteonectin immunoreactivity in the diagnosis of osteosarcoma. *Journal of Clinical Pathology*. 108: 464-473. 1997.
7. Hunzelmann N, Hafner M, Anders S, Krieg T, Nischt R. BM-40 (osteonectin, SPARC) is expressed both in the epidermal and in the dermal compartment of adult human skin. *Journal of Investigative Dermatology*. 110: 122-126. 1998.
8. Le Bail B, Faouzi S, Boussarie L, Guirouilh J, Blanc JF, Carles J, Bioulac-Sage P, Balabaud C, Rosenbaum J. Osteonectin/SPARC is overexpressed in human hepatocellular carcinoma. *Journal of Pathology*. 189: 46-52. 1999.
9. Massi D, Franchi A, Borgognoni L, Reali UM, Santucci M. Osteonectin expression correlates with clinical outcome in thin cutaneous malignant melanomas. *Human Pathology*. 30 (3): 339-344. 1999.
10. Rempel S A, Ge S and Gutiérrez J A. SPARC: a potential diagnostic marker of invasive meningiomas. *Clinical Cancer Research*. 5: 237-241. 1999.
11. Yiu GK, Chan WY, Ng SW, Chan PS, Cheung KK, Berkowitz RS, Mok SC. SPARC (secreted protein acidic and rich in cysteine) induces apoptosis in ovarian cancer cells. *American Journal of Pathology*. 159(2): 609-622. 2001.
12. Bloomston M, Ellison EC, Muscarella P, Al-Saif O, Martin EW, Melvin WS, Frankel WL. Stromal osteonectin overexpression is associated with poor outcome in patients with ampullary cancer. *Ann Surg Oncol*; 14(1):211-217. 2007.
13. Bozkurt SU, Ayan E, Bolukbasi F, Elmaci I, Pamir N, Sav A. Immunohistochemical expression of SPARC is correlated with recurrence, survival and malignant potential in meningiomas. *APMIS*; 117(9): 651-659. 2009.
14. Wang HY, Li YY, Shao Q, Hou JH, Wang F, Cai MB, Zeng YX, Shao JY. Secreted protein acidic and rich in cysteine (SPARC) is associated with nasopharyngeal carcinoma metastasis and poor prognosis. *J Transl Med*. 2012;10:27
15. Nagai MA, Gerhard R, Fregnani JH, Nonogaki S, Rierger RB, Netto MM, Soares FA. Prognostic value of NDRG1 and SPARC protein expression in breast cancer patients. *Breast Cancer Res Treat*. 2011;126:1-14
16. Witkiewicz AK, Freydin B, Chervoneva I, Potoczek M, Rizzo W, Rui H, Brody JR, Schwartz GF, Lisanti MP. Stromal CD10 and SPARC expression in ductal carcinoma in situ (DCIS) patients predicts disease recurrence. *Cancer Biol Ther*. 2010;10:391-396

17. Hsiao YH, Lien HC, Hwa HL, Kuo WH, Chang KJ, Hsieh FJ. SPARC (osteonectin) in breast tumors of different histologic types and its role in the outcome of invasive ductal carcinoma. *Breast J.* 2010;16:305-308

18. Capper D, Mittelbronn M, Goepfert B, Meyermann R, Schittenhelm J. Secreted protein, acidic and rich in cysteine (SPARC) expression in astrocytic tumour cells negatively correlates with proliferation, while vascular SPARC expression is associated with patient survival. *Neuropathol Appl Neurobiol.* 2010;36:183-197

19. Liang JF, Wang HK, Xiao H, Li N, Cheng CX, Zhao YZ, Ma YB, Gao JZ, Bai RB, Zheng HX. Relationship and prognostic significance of SPARC and VEGF protein expression in colon cancer. *J Exp Clin Cancer Res.* 2010;29:71