CD15 - Anticorpo Monoclonal anti-humano - Clone (MMA)

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1. Arber DA, Weiss LM. CD15: a review. *Appl Immunochim.* 1993, 1(1): 17-30.
2. Leong ASY, Cooper K, Leong FJWM: CD15 Manual of diagnostic antibody for Immunohistochemistry. 2th ed., London: Greenwich Medical Media, 2003. CD15, p. 83-84
3. Knowles DM. Immunophenotypic markers useful in the diagnosis and classification of hematopoietic neoplasms. In Knowles DM. *Neoplastic hematopathology.* 2th ed. Philadelphia: Lippincott, Williams & Wilkins, 2001. p. 93-226.
4. Cui L, Johkura K, Yue F, Ogiwara N, Okouchi Y, Asanuma K, Sasaki K. Spatial Distribution and Initial Changes of SSEA-1 and Other Cell Adhesion-related Molecules on Mouse Embryonic Stem Cells Before and During Differentiation. *J Histochem Cytochem.* 2004, 52 (11): 1447-1457.
5. Theilgaard-Mönch K, Raaschou-Jensen K, Schjødt K, Heilmann C, Vindeløv L, Jacobsen N, Dickmeiss E. Pluripotent and myeloid-committed CD34+ subsets in hematopoietic stem cell allografts. *Bone Marrow Transplant,* 2003, 32(12):1125-1133.
6. Sewell HF, Jaffray B, Thompson WD. Reaction of monoclonal anti Leu M1--a myelomonocytic marker (CD15)--with normal and neoplastic epithelia. *J Pathol.* 1987, 151(4):279-84.
7. McCarthy N, Simpson J, Coghill G, Kerr MA. Expression in normal adult, fetal, and neoplastic tissues of a carbohydrate differentiation antigen recognised by antigranulocyte mouse monoclonal antibodies. *J Clin Pathol.* 1985, 38(5):521-9.
8. Wasielewski R, Mengel M, Fischer R, Hansmann ML, Hübner K, Franklin J, Tesch H, Paulus U, Werner M, Diehl V, Georgii A. Classical Hodgkin's disease. Clinical impact of the immunophenotype. *Am J Pathol.* 1997, 151(4): 1123-1130.
9. Rüdiger T, Ott G, Ott MM, Müller-Deubert SM, Müller-Hermelink HK. Differential diagnosis between classic Hodgkin's lymphoma, T-cell-rich B-cell lymphoma, and paragranuloma by paraffin immunohistochemistry. *Am J Surg Pathol.* 1998, 22(10):1184-91.
10. Pileri SA, Ascani S, Leoncini L, Sabattini E, Zinzani PL, Piccaluga PP, Pileri A Jr, Giunti M, Falini B, Bolis GB, Stein H. Hodgkin's lymphoma: the pathologist's viewpoint. *J Clin Pathol.* 2002, 55(3):162-76.
11. Hall PA, D'Ardenne AJ. Value of CD15 immunostaining in diagnosing Hodgkin's disease: a review of published literature. *J Clin Pathol.* 1987, 40 (11):1298 -1304.
12. Dinand V, Malik A, Unni R, Arya LS, Pandey RM, Dawar R. Proliferative index and CD15 expression in pediatric classical Hodgkin lymphoma. *Pediatr Blood. Cancer.* 2008, 50(2):280-283.
13. Swerdlow SH, Campo E, Harris NL, Jaffe ES, Pileri SA, Stein H, Thiele J, Vardiman JW. WHO classification of tumours of haematopoietic and lymphoid tissue. 4th ed. Lyon: IARC. 2008.
14. Klebe S, Nurminen M, Leigh J, Henderson DW. Diagnosis of epithelial mesothelioma using tree-based regression analysis and a minimal panel of antibodies. *Pathology.* 2009, 41(2):140-148.
15. Yaziji H, Battifora H, Barry TS, Hwang HC, Bacchi CE, McIntosh MW, Kussick SJ, Gown AM. Evaluation of 12 antibodies for distinguishing epithelioid mesothelioma from adenocarcinoma: identification of a three-antibody immunohistochemical panel with maximal sensitivity and specificity. *Mod Pathol.* 2006, 19(4):514-523.
16. Zhou M, Roma A, Magi-Galluzzi. The usefulness of immunohistochemical markers in the differential diagnosis of renal neoplasms. *Clin Lab Med.* 2005, 25(2):247-257.
17. Pellegrini W, Bresciani G, De Zorzi A, Marocolo D, Ungari M, Facchetti. MMA monoclonal antibody is a superior anti-CD15 reagent for the diagnosis of classical Hodgkin's lymphoma? *Haematologica.* 2007, 92(5):708-709.
18. Benharroch D, Dima E, Levy A, Ohana-Malka O, Ariad S, Prinsloo I, Mejirovsky E, Sacks M, Gopas J. Differential Expression of Sialyl and Non-Sialyl-CD15 Antigens on Hodgkin-Reed-Sternberg Cells: Significance in Hodgkin's Disease. *Leuk Lymphoma.* 2000, 39 (1-2):185-194.
19. Nordic QC [en línea]: CD15. <<http://www.nordiqc.org/Epitopes/CD15/CD15.htm>> [Consulta: 8 mayo 2010]. 20. Sociedad Española de Anatomía Patológica [en línea]: Módulo de Patología Linfóide. Primera Ronda. CD15.