## CD68 - Anticorpo Monoclonal anti-humano - Clone (KP-1)

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- 1.Pulford KAF, Rigney EM, Micklem KJ, Jones M, Stross WP, Gatter KC, et al. KP1: a new monoclonal antibody that detects a monocyte/macrophage associated antigen in routinely processed tissue sections. *J Clin Pathol*; 42: 414-421. 1989.
- 2.Baldus SE; Thiele J; Park YO; Charles A; Mross C; Hanisch FG; Zirbes TK; Wickenhauser C; Fischer R. Carbohydrate and peptide antigens in macrophage populations derived from human bone marrow and milk: an immunomorphological and immunochemical analysis. *Histochemical Journal*; 27(8):630-638.1995.
- 3.Gloghini A; Rizzo A; Zanette I; Canal B; Rupolo G; Bassi P; Carbone A. KP1/CD68 expression in malignant neoplasms including lymphomas, sarcomas, and carcinomas. *American Journal of Clinical Pathology*; 103(4):425-431.1995.
- 4.Maluf HM; DeYoung BR; Swanson PE; Wick MR. Fibroma and giant cell tumor of tendon sheath: a comparative histological and immunohistological study. *Modern Pathology*; 8(2):155-159.1995.
- 5.Mazal PR;Hainfellner JA; Preiser J; Czech T; Simonitsch I; Radaszkiewicz T; Budka H. Langerhans cell histiocytosis of the hypothalamus: diagnostic value of immunohistochemistry. *Clinical Neuropathology*;15(2):87-91.1996
- 6.Menke DM; Griesser H; Araujo I; Foss HD; Herbst H; Banks PM; Stein H. Inflammatory pseudotumors of lymph node origin show macrophage-derived spindle cells and lymphocyte-derived cytokine transcripts without evidence of T-cell receptor gene rearrangements. Implications for pathogenesis and classification as an idiopathic retroperitoneal fibrosis-like sclerosing immune reaction. *American Journal of Clinical Pathology*; 105(4): 430-439.1996.
- 7.Ono T; Muso E; Suyama K; Oyama A; Matsushima H; Yashiro M; Kuwahara T; Yoshida H; Kanatsu K; Sasayama S. Intraglomerular deposition of intact cross-linked fibrin in IgA nephropathy and Henoch-Schonlein purpura nephritis. *Nephron*;74(3):522-528.1996.
- 8.Roggendorf W; Strupp S; Paulus W. Distribution and characterization of microglia/macrophages in human brain tumors. *Acta Neuropathologica*;92(3):288-293.1996.
- 9.Tetlow LC; Woolley DE. Eosinophils are an insignificant cellular component of rheumatoid synovium in patients with late stage disease: comparative distributions with mast cells and macrophages. *Annals of the Rheumatic Diseases*; 55(8):548-551.1996
- 10.Zeng L; Takeya M; Ling X; Nagasaki A; Takahashi K. Interspecies reactivities of anti-human macrophage monoclonal antibodies to various animal species. *Journal of Histochemistry and Cytochemistry*; 44(8):845-853.1996.
- 11.Zelger B; Weinlich G; Steiner H; Zelger BG; Egarter-Vigl E. Dermal and subcutaneous variants of plexiform fibrohistiocytic tumor. *American Journal of Surgical Pathology*, 1997 Feb, 21(2):235-41.
- 12.Pernick NL, DaSilva M, Gangi MD, Crissman J, Adsay V. "Histiocytic markers" in melanoma. *Mod Pathol*;12:1072-1077.1999.