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1. Lee EJ, Lee HG, Park SH, Choi EY, Park SH. CD99 type II is a determining factor for the differentiation of primitive neuroectodermal cells. *Exp Mol Med*. 2003 31;35:438-47
2. Rocchi A, Manara MC, Sciandra M, Zambelli D, Nardi F, Nicoletti G, Garofalo C, Meschini S, Astolfi A, Colombo MP, Lessnick SL, Picci P, Scotlandi K. CD99 inhibits neural differentiation of human Ewing sarcoma cells and thereby contributes to oncogenesis. *J Clin Invest*. 2010 Mar 1;120(3):668-80
3. Wick MR. Immunohistology of neuroendocrine and neuroectodermal tumors. *Semin Diagn Pathol*. 2000 Aug;17(3):194-203
4. Olsen SH, Thomas DG, Lucas DR. Cluster analysis of immunohistochemical profiles in synovial sarcoma, malignant peripheral nerve sheath tumor, and Ewing sarcoma. *Mod Pathol*. 2006 May;19(5):659-68
5. Robertson PB, Neiman RS, Worapongpaiboon S, John K, Orazi A. O13 (CD99) positivity in hematologic proliferations correlates with TdT positivity. *Mod Pathol*. 1997 Apr;10(4):277-82
6. Lucas DR, Bentley G, Dan ME, Tabaczka P, Poulik JM, Mott MP. Ewing sarcoma vs lymphoblastic lymphoma. A comparative immunohistochemical study. *Am J Clin Pathol*. 2001 Jan;115(1):11-7
7. Hamilton G, Fellinger EJ, Schratter I, Fritsch A. Characterization of a human endocrine tissue and tumor-associated Ewing's sarcoma antigen. *Cancer Res*; 48: 6127-6131. 1988.
8. Weidner N, Tjoe J. Immunohistochemical profile of monoclonal antibody O13: antibody that recognizes glycoprotein p30/32MIC2 and is useful in diagnosing Ewing's sarcoma and peripheral neuroepithelioma. *Am J Surg Pathol*; 18: 486-494. 1994.
9. Hahn JH, Kim MK, Choi EY, Kim SH, Sohn HW, Ham DI, Chung DH, Kim TJ, Lee WJ, Park CK, Ree HJ, Park SH. CD99 (MIC2) regulates the LFA-1/ICAM-1-mediated adhesion of lymphocytes, and its gene encodes both positive and negative regulators of cellular adhesion. *J Immunol*; 159: 2250-2258. 1997.
10. Wacławiczek M, Majdic O, Stulnig T, Berger M, Sunder-Plassmann R, Zlabinger GJ, Baumruker T, Stockl J, Ebner C, Knapp W, Pickl WF. CD99 engagement on human peripheral blood T cells results in TCR/CD3-dependent cellular activation and allows for Th1-restricted cytokine production. *J Immunol*; 161: 4671-4678. 1998.
11. Sung CO, Young H, Ko YH, Park S, Kim K, Kim W. Immunoreactivity of CD99 in Non-Hodgkin's Lymphoma: Unexpected Frequent Expression in ALK-positive Anaplastic Large Cell Lymphoma. *J Korean Med Sci*; 20: 952-956. 2005.
12. Buxton D, Bacchi CE, Gualco G, Weiss LM, Zuppan CW, Rowsell EH, Huang Q, Wang J. Frequent expression of CD99 in anaplastic large cell lymphoma: a clinicopathologic and immunohistochemical study of 160 cases. *Am J Clin Pathol*. 2009 Apr;131(4):574-9
13. Lee SP, Park S, Park J, Hong J, Ko YH. Clinicopathologic characteristics of CD99-positive diffuse large B-cell lymphoma. *Acta Haematol*. 2011;125(3):167-7
14. Hartel PH, Fanburg-Smith JC, Frazier AA, Galvin JR, Lichy JH, Shilo K, Franks TJ. Primary pulmonary and mediastinal synovial sarcoma: a clinicopathologic study of 60 cases and comparison with five prior series. *Mod Pathol*. 2007 Jul;20(7):760-9
15. Guo Y, Yuan F, Deng H, Wang HF, Jin XL, Xiao JC. Paranuclear dot-like immunostaining for CD99: a unique staining pattern for diagnosing solid pseudopapillary neoplasm of the pancreas. *Am J Surg Pathol*. 2011 Jun;35(6):799-806
16. Makhoul R, Schwartz AM, Williams R, Lee NH. Paranuclear dot-like immunostaining for CD99: presence in colonic adenomas and adenocarcinomas. *Am J Surg Pathol*. 2011 Nov;35(11):1749-50
17. Kim HS, Kim JS, Kim JS, Park JT, Lee MC, Juhng SW, Cho JH, Park CS. The association between CD99 and LMP-1 expression in nasopharyngeal carcinoma. *Exp Oncol*. 2006 Mar;28(1):40-3
18. Vang R, Herrmann ME, Tavassoli FA. Comparative immunohistochemical analysis of granulosa and sertoli components in ovarian sex cord-stromal tumors with mixed differentiation: potential implications for derivation of sertoli differentiation in ovarian tumors. *Int J Gynecol Pathol*. 2004 Apr;23(2):151-61
19. Milanezi F, Pereira EM, Ferreira FV, Leitão D, Schmitt FC. CD99/MIC-2 surface protein expression in breast carcinomas. *Histopathology*. 2001 Dec;39(6):578-83