Citoqueratina AE1-AE3 – Anticorpo Monoclonal anti-humano

Mouse anti-human Keratin AE1-AE3 Cocktail Antibody (Clone AE1-AE3)

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1. Moll R, Franke WW, Schiller DL, Geiger B, Krepler R: The catalog of human cytokeratins: patterns of expression in normal epithelia, tumors and cultured cells. *Cell*. 1982; 31: 11-24.
2. Tseong SC, Jarvinen MJ, Nelson WG, Huang JW, Woodcock-Mitchell J, Sun TT: Correlation of specific keratins with different types of epithelial differentiation: monoclonal antibody studies. *Cell*. 1982; 30:361-72.
3. Cooper D, Schermer A, Pruss R, Sun TT: The use of a1F, AE1, and AE3 monoclonal antibodies for the identification and classification of mammalian epithelial keratins. *Differentiation*. 1984; 28:30-5.
4. Pinkus GS, O'Connor EM, Etheridge CL, Corson JM: Optimal immunoreactivity of keratin proteins in formalin-fixed, paraffin-embedded tissue requires preliminary trypsinization. An immunoperoxidase study of various tumours using polyclonal and monoclonal antibodies. *J Histochem Cytochem*. 1985; 33:465-73.
5. Pinkus GS, Etheridge CL, O'Connor EM: Are keratin proteins a better tumor marker than epithelial membrane antigen? A comparative immunohistochemical study of various paraffin-embedded neoplasms using monoclonal and polyclonal antibodies. *Am J Clin Pathol*. 1986; 85:269-77.
6. Listrom MB, Dalton LW: Comparison of keratin monoclonal antibodies MAK-6, AE1:AE3, and CAM-5.2. *Am J Clin Pathol*. 1987 Sep;88(3):297-301.
7. Levy R, Czernobilsky B, Geiger B: Subtyping of epithelial cells of normal and metaplastic human uterine cervix, using polypeptide-specific cytokeratin antibodies. *Differentiation*. 1988; 39: 185-96.
8. Mygind H, Nielsen B, Moe D, Clausen H, Dabelsteen E, Clausen PP: Antikeratin antibodies in routine diagnostic pathology. A comparison of 10 different commercial antikeratins. *APMIS*. 1988; 96:1009-22.
9. Cosgrove M, Fitzgibbons PL, Sherrod A, Chandrasoma PT, Martin SE: Intermediate filament expression in astrocytic neoplasms. *Am J Surg Pathol*. 1989; 13:141-5.
10. Heatley MK: Cytokeratins and cytokeratin staining in diagnostic histopathology. *Histopathology*. 1996; 28: 479-83.
11. Hesse M, Magin TM, Weber K: Genes for intermediate filament proteins and the draft sequence of the human genome: novel keratin genes and a surprisingly high number of pseudogenes related to keratin genes 8 and 18. *J Cell Sci*. 2001; 114: 2569-75.
12. Porter RM, Lane EB: Phenotypes, genotypes and their contribution to understanding keratin function. *Trends Genet*. 2003; 19: 278-85.
13. Chu PG, Weiss LM: Keratin expression in human tissues and neoplasms. *Histopathology* 2002, 40, 403-439.