# Sonda Evolume FGFR3/IGH Dual Color FISH, kit

Evolume FGFR3/IGH Dual Color FISH Probe Kit

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**Nome:** Sonda Evolume FGFR3/IGH FISH Dual Color, kit

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# Sonda Evolume FGFR3/IGH Dual Color FISH, kit

Evolume FGFR3/IGH Dual Color FISH Probe Kit


**Product Name: Evolume FGFR3/IGH Dual Color FISH Probe, kit**

## Packaging Specifications

Reference	Description	Kit Size
EP-13-10971	FGFR3/IGH FISH Probe Kit	5 Tests/Kit
EP-13-10972	FGFR3/IGH FISH Probe Kit	10 Tests/Kit
EP-13-10973	FGFR3/IGH FISH Probe Kit	20 Tests/Kit

## Test Principle

### Probe Design Description

### Technique Description

Fluorescence in situ hybridization (FISH) is a molecular assay technique that uses fluorescently labeled DNA or RNA fragments or probes to hybridize via base complementary mechanism to the targeted regions of chromosomes in specimens. The probes emit fluorescence signals at a specific wavelength after excitation by a light source, and the hybridization to the chromosome can be directly observed through a microscope equipped with a proper filter set. This technique can accurately detect chromosomal abnormalities, therefore, the status of genes located at the affected regions of chromosomes can be evaluated.

## Main Components

Specification	Kit content	Main components	Volume	Quantity
5 tests/kit	FGFR3/IGH probe (FGFR3 orange probe, IGH green probe)	FGFR3/IGH probe, formamide, SSC, and dextransulfate, etc.	50 µl	1 vial
10 tests/kit			100 µl	
20 tests/kit			200 µl	
	DAPI Mounting Medium	Antifade reagent ,DAPI, and glycerin, etc.	50 µl	1 vial
			100 µl	
			200 µl	

**Notes: The components in kits of different lots should not be interchanged.**

### Material need but not provided

- Rubber cement
- Absolute ethanol
- Purified water
- 1 M HCl
- 1 N or 2 N Hidróxido de Sódio (NaOH)
- 20×SSCSolution (3 Msodium chloride, 0.3 M sodium citrate, pH 5.3)
- NP-40
- Methanol
- Glacialaceticacid
- 0.075M Potassium chloride (KCl)
- Pepsin (250 U/mg)

### Applicable instruments

- Fluorescence microscope.
- Microliter pipettor and sterile tips
- Microcentrifuge
- Water bath
- Coplin jars
- Hybridization instrument

**Notes:** Fluorescence microscope. The configuration of the required fluorescence microscope includes: 10× eye piece, and 10×, 40× and 100× objective lenses. It is recommended that before using the probe, the user should ask the filter set supplier for the details of the filter set to be used, so as to choose a filter set that is compatible with the labeled fluorescent dyes.

Orange fluorescence: excitation maximum: 552 nm, emission maximum: 576 nm.

Green fluorescence: excitation maximum: 496 nm, emission maximum: 520 nm.

DAPI: excitation maximum: 340 nm, emission maximum: 488 nm.

## Storage Conditions and Shelf Life

Storage conditions: below -15°C, sealed, and stored in the dark.

Shipping conditions: The kits must be transported under controlled temperature conditions between 2 °C and 8 °C. The shipping time must not exceed 10 days. During transport, the temperature must not exceed room temperature.

Shelf life: 12 months.

**Note:** see the production date and shelf life on the outer package

## Sample Requirements

1. Sample type: peripheral blood, or bone marrow.
2. Sample collection method: 2-3 ml (anticoagulated by heparin sodium). The sample size depends on the specific conditions of the donor.
3. Sample storage: use the sample within 2 hours after sampling, or store it at 2-8°C.

## Test Methods

### Preparation of Working Reagent

#### Ethanol Wash Solutions

Prepare v/v dilutions of 70% and 85% using 100% ethanol and purified water. Dilutions may be used for 1 week unless evaporation occurs or the solution becomes diluted due to excessive use. Store at room temperature in tightly capped containers when not in use.

### Pretreatment Solution

#### 2×SSC solution

- To prepare, add together:

100ml	20×SSCpH5.3
<u>900 mL</u>	Purified water
1000 mL	Final Volume

Mix thoroughly. Measure pH at room temperature using a pHmeter. Adjust pH to 7.0±0.2 with 1 N or 2N NaOH. Adjust volume to 1 liter with purified water. If the reagent is turbid or contaminated, it should be discarded immediately.

#### 0.4×SSC/0.3% NP-40

- To prepare, add together:

20 mL	20× SSC pH 5.3
877 mL	Purified water
<u>3 mL</u>	NP-40
1000 mL	Final volume

Mix thoroughly. Measure pH at room temperature using a pH meter. Adjust pH to 7.0±0.2 with 1 N or 2 N NaOH. Adjust volume to 1 liter with purified water. If the reagent is turbid or contaminated, it should be discarded immediately.

### **2×SSC/0.1%NP-40**

- To prepare, add together:

20 mL	20× SSC pH 5.3
849 mL	Purified water
1 mL	NP-40
1000 mL	Final volume

Mix thoroughly. Measure pH at room temperature using a pH meter. Adjust pH to  $7.0 \pm 0.2$  with 1 N or 2 N NaOH. Adjust volume to 1 liter with purified water. If the reagent is turbid or contaminated, it should be discarded immediately.

### **50 mg/ml Pepsin solution**

- To prepare, add together:

25 mg	Pepsin (250U/mg)
500 $\mu$ l	Purified water
1000 mL	Final volume

Mix thoroughly until it is fully dissolved. Dispense at 100  $\mu$ l per vial, and store the vials at  $-20 \pm 5^\circ\text{C}$ ; avoid repeated freezing and thawing.

### **Enzyme buffer (pH 2.0 HCl)**

- To prepare, add together:

10 mL	1 M HCl
900 mL	Purified water
1000 mL	Final Volume

Mix thoroughly. Measure pH at room temperature using a pH meter. Adjust pH to  $2.0 \pm 0.2$  with 1 M HCl. Adjust volume to 1 liter with purified water. If the reagent is turbid or contaminated, it should be discarded immediately.

### **Sample treatment**

**Note:** Fixative (3:1 methanol: glacial acetic acid) should be prepared fresh daily.

1. Take 5-10 ml cell samples, centrifuge them at 1500 rpm for 5 minutes, and remove the supernatant carefully.
2. Add 10 ml 0.075M KCl preheated to  $37^\circ\text{C}$  into the centrifuge tube, and mix the substances by pipetting gently.
3. Place the centrifuge tube in the  $37^\circ\text{C}$  water bath for 30 minutes.
4. Add 1 ml fresh fixative solution, mix the substances by pipetting gently, and pre-immobilize at room temperature for 8 minutes.
5. Mix the substances by pipetting and centrifuge the mixture at 1500 rpm for 8 minutes.
6. Remove the supernatant, add 10 ml fresh fixative solution to the precipitate, mix the substances by pipetting gently, and let the mixture stand at room temperature for 8 minutes.
7. Centrifuge the mixture at 1500 rpm for 8 minutes.

8. Repeat the above washing steps until the cell pellet is washed clean.

### Slide preparation

1. Take a clean object slide.
2. After resuspending the cells, take 3µl suspension and drop it to different positions on the object slide, avoiding overlap.
3. Dry the slide at room temperature.
4. Use a 20× objective lens to observe the cell density in three areas under a phase-contrast microscope; it is required that the cells have no overlap and the number of cells is 100~200; ① If the cells overlap, add an appropriate amount of fresh fixative solution to dilute the cell suspension, and after mixing, take another 3µl suspension for slide preparation; ② If the cell density is low, centrifuge the suspension, carefully remove an appropriate amount of supernatant, and after mixing, take another 3µl suspension for slide preparation; then dry the slide for observation.
5. Under the phase-contrast microscope, if there are too many cell fragments, it is necessary to do pretreatment and select the appropriate hybridization area.

**Note:** at least one extra slide is required for each sample, and the remaining cell suspension can be stored at 2-8°C for 1 month, so that slides can be re-prepared if necessary.

### Preparation of instruments

**Note:** 1 mg/ml Pepsin solution (Take 50 mg pepsin, add 50 ml enzyme buffer and mix well) should be prepared fresh daily

1. Turn on the thermostat water bath and set the temperature to 37°C. Pour 50 ml 0.1 mg/ml Pepsin solution into the Coplin jar, put the jar in the water bath, and preheat in to 37°C.
2. Before elution, set the temperature of the thermostat water bath to 68°C, pour the washing solution into the Coplin jar, and place the jar in the water bath for at least 30 min. Ensure that the temperature of the washing solution reaches 68±1°C before use.

### Testing process

#### Sample pretreatment

1. Immerse the prepared slide in 2×SSC solution at room temperature for 10 min.
2. Place the slide in 37±1°C 0.1mg/ml Pepsin solution to digest for 10 min.
3. Immerse the slide in 2×SSC solution at room temperature for 2min.
4. Dehydrate the slide in 70% ethanol, 85% ethanol and 100% ethanol in turn for 1 min each.
5. Dry the slide at room temperature.

#### Sample denaturation and hybridization

1. Place the moisturizing strip in the corresponding slot of the hybridization instrument after soaking it in purified water, and ensure that the surface of the hybridization instrument is clean and free of foreign matter.
2. The denaturation and hybridization program is set as follows: denaturation: 85°C, 5 min; hybridization: 42°C, 16-18 hours.

**Note:** The fast hybridization program: denaturation: 85°C, 5 min; hybridization: 42°C, 2 hours.

1. Take out the probes from below -15°C storage, swirl and centrifuge them after they return to room temperature, take 10 µl probes, and drop them onto the target area of the slide.

2. Place the cover slide, avoiding air bubbles, and seal the edge of the cover slide with Rubber Cement.
3. Put the sealed slide on the hybridization instrument and start the denaturation and hybridization program.

### Post-hybridization wash

1. Preheat 0.4× SSC/0.3% NP-40 to 68±1°C, and place 2× SSC/0.1% NP-40 at room temperature.
2. After the hybridization is completed, take out the cell slide, and remove the fluoro gel and cover slide from the cell slide.
3. Put the slide in 68±1°C 0.4× SSC/0.3% NP-40, and immerse and rinse it for 2 min.
4. Put the slide in 2× SSC/0.1% NP-40 at room temperature, and immerse it for 1 min.
5. Put the slide in 70% alcohol solution, and wash it for 1 min.
6. Take out the slide and dry it naturally in the dark.


### Counterstaining

1. Drop 10 µl counterstaining solution to the hybridization area (avoiding bubbles), place the cover slide, carefully remove the excess counterstaining solution with water absorbing paper, and counterstain the slide at -20±5°C in the dark for more than 20 min.
2. Select the appropriate filter set to observe the results under the fluorescence microscope.

### Result analysis

1. View the slide under a fluorescence microscope equipped with proper filter sets.
2. Place one drop of immersion oil on top of the coverslip in the area where the specimen is located. Find the specimen using the 10× objective lens and then examine the signals under the 100× lens.
3. Adjust the focus to find the observable area that has complete nucleus boundary, uniformed staining, and clear signals without overlapping nuclei.
4. In the selected region, find all signal points at different levels of the nuclei, and count the fused signals and the individual orange and green signals; when the orange (O) and green (G) signal points overlap or the distance between the two signal points is less than one signal point diameter, the signals are counted as a fused signal (yellow, F); the nuclei without signals or with weak signals are not counted.

## Interpretation of common signal type

FISH scenario 2O2G (2 orange / 2 green signals)	
FISH scenario 2F1O1G (2 yellow / 1 orange / 1 green signals)	

### Threshold Determination for FISH Negative Specimen

Twenty specimens without the known genetic abnormalities to the FISH probe set were randomly selected to prepare control slides. After hybridization, count 200 cells from each control slide and record any types of signal patterns that are considered to be FISH positive. Calculate the percentage of positive signals observed and determine the standard deviation. The threshold for FISH negative specimen is set as the average value of percentage plus 3 times of the standard deviation. The control or FISH negative specimens are usually chosen to be the same type of tissue or cell as that of the intended target of detection by FISH.

### Importance of Threshold Determination

The threshold must be set when probes are used for the first time since this value will be served as the reference to distinguish FISH positive and negative specimens. If the experimental procedures are altered such as the method of specimen pre-treatment or change of equipment, the threshold may vary therefore must be reset under the newly established experimental conditions.

## Interpretation of Test Results

### Common problems and solutions in the testing process:

Problem	Possible cause	Recommended solution
No signal or weak signal	Insufficient denaturation of samples and probes	Ensure that the temperature of the slide is 85±1°C during the denaturation; Extend the slide denaturation time by 2-4 min.
	No probe added	Fully thaw the probes, and ensure that the probe reagent is drawn by the pipette.
	Insufficient quantity of probes	Ensure that the probes reach room temperature before use, and ensure that the probe reagent is drawn by the pipette.
	Insufficient drying of slide	Before adding the probes onto the slide, ensure that the ethanol solution on the slide has completely volatilized.
	Too fast drying of probes	Cover the target area by the cover slide immediately after the probes are added; for the elution, remove the cover slide on only one slide at a time, and immerse the slide in the washing solution immediately before the cover slide on the next slide is removed.
	Bolhas sob a lâmina durante a hibridização	Cover the surface of the probes with the cover slide, and gently squeeze it to eliminate the air bubbles.
	Inappropriate hybridization conditions	Ensure the compliance with the specified hybridization time and temperature; do not leave any gap when sealing the slide with fluoro gel; Adjust the hybridization time and humidity as appropriate.
	Incorrect washing solution or elution conditions	Ensure that the washing solution is prepared in accordance with the IFU; ensure that the temperature of the washing solution reaches the temperature specified in the elution step; remove the cover slide before immersing the slide in the washing solution.
	Improper storage of probes or sample slides	Store the probes below -15°C in the dark; dry the unhybridized slides for long-term storage at -20±5°C or short-term storage (generally no more than two weeks) at room temperature; Store the hybridized slides at -20±5°C in the dark.
	Inappropriate selection of filter set for observation	Use the correct filter set to observe the probe fluorescence.
Inappropriate microscope structure and objective lens for observing FISH samples, or damage of the filter set	Please contact the microscope manufacturer.	
Excessively intense background of slides	Insufficient washing of slides before sample preparation	Immerse the slide in absolute ethanol and wipe it dry with water absorbing paper before dropping the reagent.
	Inadequate elution after the hybridization	Ensure that the washing solution is prepared correct in accordance with the IFU; operate in accordance with the IFU; ensure that the temperature of the washing solution is correct; remove the cover slide and repeat the elution step.
	Long-term use or improper storage of washing solution	Ensure that the washing solution is stored at 2-25°C, and discard the reagent if it is turbid or contaminated.
	Too high or too low hybridization humidity	Adjust the hybridization humidity to the optimum.
Too weak counterstaining	Too weak counterstaining	Remove the cover slide and immerse the slide in the washing solution for 5 min at room temperature. Place the slide in 70%, 85% and 100% ethanol solutions for 1 min each for gradient dehydration, and then perform counterstaining.
	Expiration or excessive light exposure of staining solution	Ensure that the staining solution is stored below -15°C in the dark; ensure that the staining solution is not expired.


## Limitations of Test Methods

1. This kit is based on the fluorescence in situ hybridization technique for the FGFR3/IGH gene Dual Fusion Translocation.
2. This kit is applicable to peripheral blood, or bone marrow.

## Precautions

1. This kit is only for research.
2. During the operation of this kit, it is necessary to wear latex gloves to avoid the contact of the reagent with the skin. In the case of accidental contact, rinse immediately with plenty of water.
3. Samples should not be exposed to acid and alkali or be overheated, which may damage DNA and lead to failure of FISH test.
4. All components in the kit should be used within the shelf life.
5. The staining solution contains DAPI, which is a mutagen. Inhalation, ingestion or contact with the skin should be avoided.
6. The probe contains a teratogen called formamide. Contact of it with the skin and mucous membranes should be avoided.
7. In order to obtain ideal results, it is necessary to ensure that the reagents are properly prepared and stored in accordance with the IFU.

## Explanation of Signs

	Upper limit of temperature		Keep Dry
	Do Not Use if Package is Damage		Consult instructions for use
	Do not reuse		Keep away from sunlight
	Use By		Lot Number
	Reference Number		Manufacturer
	Date of Manufacture		Contains sufficient for <n> tests