# Sonda Evolume SS18(SYT) Break-Apart FISH, kit

Evolume SS18 (SYT) Break-Apart FISH Probe Kit

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**Nome:** Sonda Evolume SS18 (SYT) Break Apart FISH, kit

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Referência	Descrição	Tamanho Kit
EP-13-11231	Sonda SS18 (SYT) Break Apart FISH, kit	5 Testes/Kit
EP-13-11232	Sonda SS18 (SYT) Break Apart FISH, kit	10 Testes/Kit
EP-13-11233	Sonda SS18 (SYT) Break Apart FISH, kit	20 Testes/Kit

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Evolume SS18 (SYT) Break-Apart FISH Probe Kit

For Research Use Only (RUO)

## Product Name: Evolume SS18 (SYT) Break-Apart FISH Probe Kit

### Packaging Specifications


Reference	Description	Kit Size
EP-13-11231	SS18(SYT) BreakApart FISH Probe Kit	5 Tests/Kit
EP-13-11232	SS18(SYT) BreakApart FISH Probe Kit	10 Tests/Kit
EP-13-11233	SS18(SYT) BreakApart FISH Probe Kit	20 Tests/Kit

### Intended Use

The SS18(SYT) BreakApart FISH Probe Kit is a test to detect rearrangements involving the SS18 gene via fluorescence in situ hybridization (FISH) in specimen to provide the physician with supporting information for evaluation.

### Test Principle

#### Probe Design Description

### Technique Description

Fluorescence in situ hybridization (FISH) is a molecular assay technique that uses fluorescently labeled DNA or RNA fragments or probes to hybridize via base complementary mechanism to the targeted regions of chromosomes in specimens. The probes emit fluorescence signals at a specific wavelength after excitation by a light source, and the hybridization to the chromosome can be directly observed through a microscope equipped with a proper filter set. This technique can accurately detect chromosomal abnormalities, therefore, the status of genes located at the affected regions of chromosomes can be evaluated.

## Main Components

Specification	Kit content	Main components	Volume	Quantity
5 tests/kit	SS18 probe (SS18,5' orange probe, SS18,3' green probe)	SS18 probe formamide, SSC, and dextran sulfate, etc.	50 µl	1 vial
10 tests/kit			100 µl	
20 tests/kit			200 µl	
	DAPI Mounting Medium	Antifade reagent ,DAPI, and glycerin,etc.	50 µl	1 vial
			100 µl	
			200 µl	

**Notes: 1.The components in kits of different lots should not be interchanged.**

### 2.Experimental reagents not provided in the kit but necessary:

- Rubber cement
- Deparaffinization solution
- Absolute ethanol
- Purified water
- 0.5M EDTA pH8.0
- 2M Tris-HCL pH8.0
- 1 M HCl
- 1 N or 2 N Hidróxido de Sódio (NaOH)
- 20×SSC Solution (3 M sodium chloride, 0.3 M sodium citrate, pH 5.3)
- NP-40
- Sodium thiocyanate
- Pepsin (250 U/mg)

### Storage Conditions and Shelf Life

Storage conditions: below -15°C, sealed, and stored in the dark.

Shipping conditions: The kits must be transported under controlled temperature conditions between 2 °C and 8 °C. The shipping time must not exceed 10 days. During transport, the temperature must not exceed room temperature.

Shelf life: 12 months.

**Note:** see the production date and shelf life on the outer package

### Applicable Instruments

- Fluorescence microscope.
- Microliter pipettor and sterile tips
- Microcentrifuge
- Waterbath
- Coplinjars

- Hybridization instrument

**Notes:** Fluorescence microscope. The configuration of the required fluorescence microscope includes: 10× eyepiece, and 10×, 40× and 100× objective lenses. It is recommended that before using the probe, the user should ask the filter set supplier for the details of the filter set to be used, so as to choose a filter set that is compatible with the labeled fluorescent dyes.

Orange fluorescence: excitation maximum: 552 nm, emission maximum: 576 nm.

Green fluorescence: excitation maximum: 496 nm, emission maximum: 520 nm.

DAPI: excitation maximum: 340 nm, emission maximum: 488 nm.

## Sample Requirements

1. Sample type: formalin-fixed, paraffin-embedded tissue sections, it is recommended that paraffin samples less than two years be selected to ensure the accuracy of the test.
2. Sample collection method: the fresh tissues should be fixed in 10% neutral formalin buffer for 24-48 hours within 0.5-1 hour after being removed. The tissue size should not exceed 0.5 cm<sup>3</sup> and the recommended section thickness is 3 µm- 5 µm.
3. Sample storage: the prepared paraffin sections are recommended to be used immediately, or stored at -20±5°C.

## Test Methods

### Preparation of Working Reagen

#### Ethanol Wash Solutions

Prepare v/v dilutions of 70% and 85% using 100% ethanol and purified water. Dilutions may be used for 1 week unless evaporation occurs or the solution becomes diluted due to excessive use. Store at room temperature in tightly capped containers when not in use.

#### Pretreatment Solution

- To prepare, add together:
 

80 g	Sodium thiocyanate
800 mL	Purified water
1000 mL	Final Volume

Mix thoroughly. Adjust volume to 1 liter with purified water. Filter through 0.45 µm pore filtration unit. Discard used solution at the end of each day. If the reagent is turbid or contaminated, it should be discarded immediately.

#### 1× EDTA Solution (10 mM EDTA, 10 mM Tris-HCL)

- To prepare, add together:
 

20 mL	0.5 M EDTA pH8.0
5 mL	2 M Tris-HCL pH 8.0
975 mL	Purified water
1000 mL	Final Volume

## 2× SSC /0.3%NP-40 Washing Solution

- To prepare, add together:

100 mL	20× SSC pH 5.3
847 mL	Purified water
<u>3 mL</u>	NP-40
1000 mL	Final volume

Mix thoroughly. Measure pH at room temperature using a pH meter. Adjust pH to  $7.0 \pm 0.2$  with 1 N or 2 N NaOH. Adjust volume to 1 liter with purified water. If the reagent is turbid or contaminated, it should be discarded immediately.

## Enzyme buffer (pH 2.0 HCl)

- To prepare, add together:

10 mL	1 M HCl
900 mL	Purified water
1000 mL	Final Volume

Mix thoroughly. Measure pH at room temperature using a pH meter. Adjust pH to  $2.0 \pm 0.2$  with 1 M HCl. Adjust volume to 1 liter with purified water. If the reagent is turbid or contaminated, it should be discarded immediately.

## Preparation of instruments

**Note:** 2 mg/ml Pepsin solution (Take 100 mg pepsin, add 50 ml enzyme buffer and mix well) should be prepared fresh daily

- Turn on the thermostat water bath and set the temperature to  $90^{\circ}\text{C}$ . Pour the Pretreatment Solution into the Coplin jar, put the jar in the water bath, and preheat it to  $90^{\circ}\text{C}$ .
- Turn on the thermostat water bath and set the temperature to  $37^{\circ}\text{C}$ . Pour 50 ml 2 mg/ml Pepsin solution into the Coplin jar, put the jar in the water bath, and preheat it to  $37^{\circ}\text{C}$ .
- Before post-hybridization wash, set the temperature of the thermostat water bath to  $68^{\circ}\text{C}$ , pour the 2× SSC/0.3% NP-40 Washing Solution into the Coplin jar, and place the jar in the water bath for at least 30 minutes. Ensure that the temperature of the washing solution reaches  $68 \pm 1^{\circ}\text{C}$  before use. Pour the 2× SSC/0.3% NP-40 Washing Solution into another Coplin jar and keep room temperature.

## Testing process

### Sample pretreatment

- Preheat the slide at  $65^{\circ}\text{C}$  for 5 minutes in an oven or using the hybridization instrument or hybridizer, recommended ARA (EP-31-20147/EasyPath) or equivalent.
- Immediately insert the preheated slide into a slide holder immersed in the deparaffinization solution (isoparaffin type) or xylene and incubate at room temperature for 15 minutes.
- After incubation, transfer the slide holder into freshly prepared deparaffinization solution and incubate for additional 15 minutes.
- After deparaffinization, remove the slide holder and place it on paper towel for 2 minutes to drain extra liquid.
- Dehydrate the slide in absolute ethanol twice for 5 minutes each.

6. Rehydrate the slide in 85% ethanol, 70% ethanol and purified water in sequence for 1 minute each
7. Take out the slide and use the water absorbing paper to remove the residual solution from the slide.
8. Place the slide in the pretreatment solution preheated to 90°C and incubate for 20 minutes.
9. Take out the slide and immerse it in purified water for 1 minute. Then, take it out and remove the residual liquid.
10. Immerse the slide in 1 mg/ml Pepsin solution preheated to 37°C and incubate for 25 minutes.
11. After the enzymatic digestion, immerse the slide in purified water for 1 minute.
12. Dehydrate the slide in 70% ethanol, 85% ethanol and 100% ethanol in sequence for 1 minute each.
13. Dry the slide at room temperature.

### Sample denaturation and hybridization

1. Warm the probes to room temperature, swirl and mix well, briefly spin the probes for 2- 3 seconds in a microcentrifuge, and drop 10 µl probes on top of the specimen. specimen. Place the coverslip (18mm x 18mm), avoid air bubbles, and seal the slide with Rubber cement.
2. Place the slide into a hybridizer and run the hybridization program as the following: denaturation at 85°C for 5 minutes and hybridization at 42°C for 16-18 hours.

**Note:** The fast hybridization program: denaturation: 85°C, 5 min; hybridization: 42°C, 2 hours.

### Post-hybridization wash

1. Carefully remove the Rubber Cement.
2. Place the slide into 2× SSC/0.3% NP-40 Washing Solution in the Coplin jar at room temperature to wash off the coverslip (2- 5 minutes).
3. Take out the slide and rinse it for 2 minutes with the 2× SSC/0.3% NP-40 Washing Solution preheated to 68±1°C.
4. Take out the slide, absorb the residual solution, put the slide in 70% ethanol at room temperature, and rinse it for 1 minute.
5. Air-dry the slide naturally in the dark for use later.

### Counterstaining

1. Drop 10 µl counterstaining solution to the hybridization area (avoiding bubbles), place the coverslip, carefully remove the excess counterstaining solution with absorbing paper, and counterstain the slide in the dark for 10 minutes.


**Note:** The counterstained slide should be examined immediately using a fluorescence microscope or stored at -15°C or lower in the dark for later observation. The stored slide should be warmed to room temperature prior to viewing under the microscope.

## Slide viewing and data collection

1. View the slide under a fluorescence microscope equipped with proper filter sets.
2. Place one drop of immersion oil on top of the coverslip in the area where the specimen is located. Find the specimen using the 10× objective lens and then examine the signals under the 100× lens.
3. Adjust the focus to find the observable area that has complete nucleus boundary, uniformed staining, and clear signals without overlapping nuclei
4. In the selected region, find all signal points at different levels of the nuclei, and count the fused signals and the individual orange and green signals; when the orange (O) and green (G) signal points overlap or the distance between the two signal points is less than one signal point diameter, the signals are counted as a fused signal (yellow, F); the nuclei without signals or with weak signals are not counted

## Interpretation of common signal types

## Threshold Determination for FISH Negative Specimen

Twenty specimens without the known genetic abnormalities to the FISH probe set were randomly selected to prepare control slides. After hybridization, count 200 cells from each control slide and record any types of signal patterns that are considered to be FISH positive. Calculate the percentage of positive signals observed and determine the standard deviation. The threshold for FISH negative specimen is set as the average value of percentage plus 3 times of the standard deviation. The control or FISH negative specimens are usually chosen to be the same type of tissue or cell as that of the intended target of detection by FISH.

## Importance of Threshold Determination

The threshold must be set when probes are used for the first time since this value will be served as the reference to distinguish FISH positive and negative specimens. If the experimental procedures are altered such as the method of specimen pre-treatment or change of equipment, the threshold may vary therefore must be reset under the newly established experimental conditions.

## Interpretation of Test Results

### Common problems and solutions in the testing process:

Problem	Possible cause	Recommended solution
No signal or weak signal	Insufficient denaturation of samples and probes	Ensure that the temperature of the slide is $85 \pm 1^\circ\text{C}$ during the denaturation; Extend the slide denaturation time by 2-4 min.
	No probe added	Fully thaw the probes, and ensure that the probe reagent is drawn by the pipette.
	Insufficient quantity of probes	Ensure that the probes reach room temperature before use, and ensure that the probe reagent is drawn by the pipette.
	Insufficient drying of slide	Before adding the probes onto the slide, ensure that the ethanol solution on the slide has completely volatilized.
	Too fast drying of probes	Cover the target area by the cover slide immediately after the probes are added; for the elution, remove the cover slide on only one slide at a time, and immerse the slide in the washing solution immediately before the cover slide on the next slide is removed.
	Bolhas sob a lâmina durante a hibridização	Cover the surface of the probes with the cover slide, and gently squeeze it to eliminate the air bubbles.
	Inappropriate hybridization conditions	Ensure the compliance with the specified hybridization time and temperature; do not leave any gap when sealing the slide with fluoro gel; Adjust the hybridization time and humidity as appropriate.
	Incorrect washing solution or elution conditions	Ensure that the washing solution is prepared in accordance with the IFU; ensure that the temperature of the washing solution reaches the temperature specified in the elution step; remove the cover slide before immersing the slide in the washing solution.
	Improper storage of probes or sample slides	Store the probes below $-15^\circ\text{C}$ in the dark; dry the unhybridized slides for long-term storage at $-20 \pm 5^\circ\text{C}$ or short-term storage (generally no more than two weeks) at room temperature; Store the hybridized slides at $-20 \pm 5^\circ\text{C}$ in the dark.
	Inappropriate selection of filter set for observation	Use the correct filter set to observe the probe fluorescence.
Excessively intense background of slides	Inappropriate microscope structure and objective lens for observing FISH samples, or damage of the filter set	Please contact the microscope manufacturer.
	Insufficient washing of slides before sample preparation	Immerse the slide in absolute ethanol and wipe it dry with water absorbing paper before dropping the reagent.
	Inadequate elution after the hybridization	Ensure that the washing solution is prepared correct in accordance with the IFU; operate in accordance with the IFU; ensure that the temperature of the washing solution is correct; remove the cover slide and repeat the elution step.
	Long-term use or improper storage of washing solution	Ensure that the washing solution is stored at $2-25^\circ\text{C}$ , and discard the reagent if it is turbid or contaminated.
Too weak counterstaining	Too high or too low hybridization humidity	Adjust the hybridization humidity to the optimum.
	Too weak counterstaining	Remove the cover slide and immerse the slide in the washing solution for 5 min at room temperature. Place the slide in 70%, 85% and 100% ethanol solutions for 1 min each for gradient dehydration, and then perform counterstaining.
	Expiration or excessive light exposure of staining solution	Ensure that the staining solution is stored below $-15^\circ\text{C}$ in the dark; ensure that the staining solution is not expired.


## Limitations of Test Methods

1. This kit is based on the fluorescence in situ hybridization technique for the SS18 gene break.
2. This reagent is applicable to paraffin tissue sections, and adjust the pretreatment process according to the thickness of the tissue sections or the type of samples.

## Precautions

1. This kit is only for research.
2. During the operation of this kit, it is necessary to wear latex gloves to avoid the contact of the reagent with the skin. In the case of accidental contact, rinse immediately with plenty of water.
3. Samples should not be exposed to acid and alkali or be overheated, which may damage DNA and lead to failure of FISH test.
4. All components in the kit should be used within the shelf life.
5. The staining solution contains DAPI, which is a mutagen. Inhalation, ingestion or contact with the skin should be avoided.
6. The probe contains a teratogen called formamide. Contact of it with the skin and mucous membranes should be avoided.
7. In order to obtain ideal results, it is necessary to ensure that the reagents are properly prepared and stored in accordance with the IFU.

## Explanation of Signs

	Upper limit of temperature		Keep Dry
	Do Not Use if Package is Damage		Consult instructions for use
	Do not reuse		Keep away from sunlight
	Use By		Lot Number
	Reference Number		Manufacturer
	Date of Manufacture		Contains sufficient for <n> tests