Sonda MYC FISH, kit

Probe MYC FISH, kit

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Nome: Sonda MYC FISH, kit

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
INTENDED USE

The MYC FISH Probe Kit is designed to detect amplification of the MYC gene via fluorescence in situ hybridization (FISH) in formalin-fixed paraffin-embedded (FFPE) non small cell lung cancer (NSCLC) tissue specimens, to provide the physician with supporting information for diagnosis.

Test Principle

Probe Design Description

The MYC/CEN8 probe is located on chromosome 8. A 1200Kb probe covering the MYC gene region is labeled with an orange dye, and part of chromosome 8 (CEN8) is labeled with a green dye. The MYC gene marker region is located on 8q24.21, which has high specificity and will not hybridize with other chromosome centromeres to produce hybrid spots.

Technique Description

Fluorescence in situ hybridization (FISH) is a molecular assay technique that uses fluorescently labeled DNA or RNA fragments or probes to hybridize via base complementary mechanism to the targeted regions of chromosomes in specimens. The probes emit fluorescence signals at a specific wavelength after excitation by a light source, and the hybridization to the chromosome can be directly observed through a microscope equipped with a proper filter set. This technique can accurately detect chromosomal abnormalities, therefore, the status of genes located at the affected regions of chromosomes can be evaluated.

Main Components

| Specification | Kit content | Main components | Volume | Quantity |
|---------------|-----------------------------------------------------------|---------------------------------------------------------------------|--------|----------|
| 5 tests/kit | MYC/CEN8 probe (MYC orange probe, CEN8 green probe) | MYC/CEN8probe, COTIDNA, formamide, SSC, and dextransulfate, etc. | 50 µl | 1 vial |
| 10 tests/kit | | | 100 µl | |
| 20 tests/kit | DAPI Mounting Medium | Antifade reagent ,DAPI, and glycerin,etc. | 200 µl | |
| | | | 50 µl | 1 vial |
| | | | 100 µl | |
| | | | 200 µl | |

Notes: The components in kits of different lots should not be interchanged.

Material need but not provided

- Rubber cement
- Deparaffinization solution
- Absolute ethanol
- Purified water
- 1 M HCl
- 1 N or 2 N Hidróxido de Sódio (NaOH)
- 20×SSCSolution (3 M sodium chloride, 0.3 M sodium citrate, pH 5.3)
- NP-40
- Sodiumthiocyanate
- Pepsin (250 U/mg)

Applicable instruments

- Fluorescence microscope.
- Microliter pipettor and sterile tips
- Microcentrifuge
- Water bath
- Coplin jars
- Hybridization instrument

Notes: Fluorescence microscope. The configuration of the required fluorescence microscope includes: 10× eye piece, and 10×, 40× and 100× objective lenses. It is recommended that before using the probe, the user should ask the filter set supplier for the details of the filter set to be used, so as to choose a filter set that is compatible with the labeled fluorescent dyes.

Orange fluorescence: excitation maximum: 552 nm, emission maximum: 576 nm.

Green fluorescence: excitation maximum: 496 nm, emission maximum: 520 nm.

DAPI: excitation maximum: 340 nm, emission maximum: 488 nm.

Storage Conditions and Shelf Life

Storage conditions: below-15°C, sealed, and stored in the dark.

Shipping conditions: The kits must be transported under controlled temperature conditions between 2 °C and 8 °C. The shipping time must not exceed 10 days. During transport, the temperature must not exceed room temperature.

Shelf life: 12 months.

Note: see the production date and shelf life on the outer package

Sample Requirements

1. Sample type: formalin-fixed,paraffin-embedded tissue sections, it is recommended that paraffin samples less than two years be selected to ensure the accuracy of the test.
2. Sample collection method: the fresh tissues should be fixed in 10% neutral formalin buffer for 24-48 hours within 0.5-1 hour after being removed. The tissue size should not exceed 0.5 cm³ and the recommended section thickness is 3 µm- 5 µm.
3. Sample storage: the prepared paraffin sections are recommended to be used immediately, or stored at -20±5°C.

Test Methods

A - Preparation of Working Reagen

Ethanol Wash Solutions

Prepare v/v dilutions of 70% and 85% using 100% ethanol and purified water. Dilutions may be used for 1 week unless evaporation occurs or the solution becomes diluted due to excessive use. Store at room temperature in tightly capped containers when not in use.

Pretreatment Solution

- To prepare, add together:

| | |
|---------|--------------------|
| 80 g | Sodium thiocyanate |
| 800 mL | Purified water |
| 1000 mL | FinalVolume |

Mix thoroughly. Adjust volume to 1 liter with purified water. Filter through 0.45 µm pore filtration unit. Discard used solution at the end of each day. If the reagent is turbid or contaminated, it should be discarded immediately.

2× SSC /0.3%NP-40 Washing Solution

- To prepare, add together:

| | |
|---------|----------------|
| 100 mL | 20× SSC pH 5.3 |
| 847 mL | Purified water |
| 3 mL | NP-40 |
| 1000 mL | Final volume |

Mix thoroughly. Measure pH at room temperature using a pH meter. Adjust pH to 7.0±0.2 with 1 N or 2 N NaOH. Adjust volume to 1 liter with purified water. If the reagent is turbid or contaminated, it should be discarded immediately.

Enzyme buffer (pH 2.0 HCl)

To prepare, add together:

| | |
|--------------|----------------|
| 10 mL | 1 M HCl |
| <u>900mL</u> | Purified water |
| 1000 mL | Final Volume |

Mix thoroughly. Measure pH at room temperature using a pH meter. Adjust pH to 2.0 ± 0.2 with 1 M HCl. Adjust volume to 1 liter with purified water. If the reagent is turbid or contaminated, it should be discarded immediately.

B - Preparation of instruments

Note: 2 mg/ml Pepsin solution (Take 100 mg pepsin, add 50ml enzyme buffer and mixwell) should be prepared fresh daily.

1. Turn on the thermostat water bath and set the temperature to 80°C. Pour the Pretreatment Solution into the Coplin jar, put the jar in the water bath, and preheat it to 80°C.
2. Turn on the thermostat water bath and set the temperature to 37°C. Pour 50 ml 1 mg/ml Pepsin solution into the Coplin jar, put the jar in the water bath, and preheat it to 37°C.
3. Before post-hybridization wash, set the temperature of the thermostat water bath to 68°C, pour the 2× SSC/0.3% NP-40 Washing Solution into the Coplin jar, and place the jar in the water bath for at least 30 minutes. Ensure that the temperature of the washing solution reaches $68 \pm 1^\circ\text{C}$ before use. Pour the 2× SSC/0.3% NP-40 Washing Solution into another Coplin jar and keep room temperature.

Testing process

Sample pretreatment

1. Preheat the slide at 65°C for 5 min on the hybridization instrument (Recommended Hybridizer: Ara EP-31-20147/EasyPath), and deparaffinize with the transparent deparaffinization solution (isoparaffintype) or xylene twice, 15 min each time.
2. Dehydrate the slide in absolute ethanol twice, 5 min each time.
3. Rehydrate the slide in 85% ethanol, 70% ethanol and purified water in gradients for 1 min each.
4. Take out the slide and use the water absorbing paper to remove the residual solution from the slide.
5. Place the slide in the pretreatment solution preheated to 80°C and maintain for 30 min.
6. Take out the slide and immerse it in purified water for 1 min; then, take it out and remove the residual solution from it.
7. Immerse the slide in 2 mg/ml Pepsin solution preheated to 37°C for 25 min.
8. After digestion, immerse the slide in purified water for 1 min.
9. Dehydrate the slide in 70% ethanol, 85% ethanol and 100% ethanol in turn for 1 min each.
10. Dry the slide at room temperature for use later.

Sample denaturation and hybridization

1. Preheat the probes to the room temperature, swirl and mix them well, centrifuge the mixture for 2s-3s in a microcentrifuge, and drop 10 µl probes to the hybridization area.
2. Place the cover slide (18mmx18mm), avoiding air bubbles, and seal the slide with rubber cement.

3. Set the hybridization instrument program: denaturation: 85°C, 5 min; hybridization: 42°C, 16-18 hours.

Note: The fast hybridization program: denaturation: 85°C, 5 min; hybridization: 42°C, 2 hours.

Post-hybridization wash

1. Carefully remove the Rubber Cement.
2. Place the slide in the room temperature 2×SSC/0.3%NP-40 Washing Solution in the Coplin jar to wash off the cover slide (2min-5min).
3. Take out the slide and rinse it for 2 min with the 2×SSC/0.3%NP-40 Washing Solution preheated to 68±1°C.
4. Take out the slide and rinse it for 1 min with the 2×SSC/0.3%NP-40 Washing Solution preheated at the room temperature
5. Take out the slide, absorb the residual solution, put the slide in 70% ethanol at room temperature, and rinse it for 1 min.
6. Air-dry the slide naturally in the dark for use later.


Counterstaining

1. Drop 10µl counterstaining solution to the hybridization area (avoiding bubbles), place the cover slide, carefully remove the excess counterstaining solution with water absorbing paper, and counter stain the slide at -20±5°C in the dark for more than 20 min.
2. Select the appropriate filter set to observe the results under the fluorescence microscope

Slide viewing and data collection

1. Select the appropriate filter set to observe the results under the fluorescence microscope.
2. Find the area under a 10× objective lens and count the number under a 100× objective lens.
3. Adjust the focal length to find the suitable area, requiring complete nucleus boundary, uniform staining by the counterstaining solution, no overlap of nuclei, and clear signal points.
4. In the selected region, find all signal points at different levels of the nuclei, and count the individual orange and green signals; the nuclei without signals or with weak signals are not counted.

Interpretation of common signal type

Threshold determination for FISH negative specimen

Twenty specimens without the known genetic abnormalities to the FISH probe set are required to be randomly selected to prepare control slides. After hybridization, count 200 cells from each control slide and record any types of signal patterns that are considered to be FISH positive. Calculate the percentage of positive signals observed and determine the standard deviation. The threshold for FISH-negative specimens is set as the average value of percentage plus 3 times the standard deviation. The control or FISH-negative specimens are usually chosen to be the same type of tissue or cell as that of the intended target of detection by FISH.

Importance of threshold determination

The threshold must be set when probes are used for the first time since this value will serve as the reference to distinguish FISH-positive and negative specimens. If the experimental procedures are altered, such as the method of specimen pre-treatment or a change of equipment, the threshold may vary; therefore must be reset under the newly established experimental conditions.

Interpretation of Test Results

Common problems and solutions in the testing process:

| Problem | Possible cause | Recommended solution |
|---------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Nosignal or weak signal | Insufficient denaturation of samples and probes | Ensure that the temperature of the slide is $85\pm 1^{\circ}\text{C}$ during the denaturation; Extend the slide denaturation time by 2-4 min. |
| | No probe added | Fully thaw the probes, and ensure that the probe reagent is drawn by the pipette. |
| | Insufficient quantity of probes | Ensure that the probes reach room temperature before use, and ensure that the probe reagent is drawn by the pipette. |
| | Insufficient drying of slide | Before adding the probes onto the slide, ensure that the ethanol solution on the slide has completely volatilized. |
| | Too fast drying of probes | Cover the target area by the cover slide immediately after the probes are added; for the elution, remove the cover slide on only one slide at a time, and immerse the slide in the washing solution immediately before the cover slide on the next slide is removed. |
| | Bubbles under the coverslip during hybridization | Cover the surface of the probes with the cover slide, and gently squeeze it to eliminate the air bubbles. |
| | Inappropriate hybridization conditions | Ensure the compliance with the specified hybridization time and temperature; do not leave any gap when sealing the slide with fluoro gel; Adjust the hybridization time and humidity as appropriate. |
| | Incorrect washing solution or elution conditions | Ensure that the washing solution is prepared in accordance with the IFU; ensure that the temperature of the washing solution reaches the temperature specified in the elution step; remove the cover slide before immersing the slide in the washing solution. |
| | Improper storage of probes or sample slides | Store the probes below -15°C in the dark; dry the unhybridized slides for long-term storage at $-20\pm 5^{\circ}\text{C}$ or short-term storage (generally no more than two weeks) at room temperature; Store the hybridized slides at $-20\pm 5^{\circ}\text{C}$ in the dark. |
| | Inappropriate selection of filter set for observation | Use the correct filter set to observe the probe fluorescence. |
| Inappropriate microscope structure and objective lens for observing FISH samples, or damage of the filter set | Please contact the microscope manufacturer. | |
| Excessively intense background of slides | Insufficient washing of slides before sample preparation | Immerse the slide in absolute ethanol and wipe it dry with water absorbing paper before dropping the reagent. |
| | Inadequate elution after the hybridization | Ensure that the washing solution is prepared correct in accordance with the IFU; operate in accordance with the IFU; ensure that the temperature of the washing solution is correct; remove the cover slide and repeat the elution step. |
| | Long-term use or improper storage of washing solution | Ensure that the washing solution is stored at $2-25^{\circ}\text{C}$, and discard the reagent if it is turbid or contaminated. |
| | Too high or too low hybridization humidity | Adjust the hybridization humidity to the optimum. |
| Too weak counterstaining | Too weak counterstaining | Remove the cover slide and immerse the slide in the washing solution for 5 min at room temperature. Place the slide in 70%, 85% and 100% ethanol solutions for 1 min each for gradient dehydration, and then perform counterstaining. |
| | Expiration or excessive light exposure of staining solution | Ensure that the staining solution is stored below -15°C in the dark; ensure that the staining solution is not expired. |


Limitations of Test Methods

1. This kit is based on the fluorescence in situ hybridization technique for the MYC gene amplification.
2. This reagent is applicable to paraffin tissue sections, and adjust the pretreatment process according to the thickness of the tissue sections or the type of samples.

Precautions

1. This kit is only for research.
2. Do not reuse the kit components.
3. The clinical significance of the test result of this kit should be evaluated in conjunction with the patient's past medical history and do their examination results.
4. During the operation of this kit, it is necessary to wear latex gloves to avoid the contact of the reagent with the skin. In the case of accidental contact, rinse immediately with plenty of water.
5. Samples should not be exposed to acid and alkali or be overheated, which may damage DNA and lead to failure of FISH test.
6. All components in the kit should be used within the shelf life.
7. The staining solution contains DAPI, which is a mutagen. Inhalation, ingestion or contact with the skin should be avoided.
8. The probe contains a teratogen called formamide. Contact of it with the skin and mucous membranes should be avoided.
9. In order to obtain ideal results, it is necessary to ensure that the reagents are properly prepared and stored in accordance with the IFU.
10. The discarded samples and wastes in the testing process should be collected and disposed of as medical waste.

Explanation of Signs

| | | | |
|-------------------------------------------------------------------------------------|---------------------------------|-------------------------------------------------------------------------------------|-----------------------------------|
|  | Upper limit of temperature |  | Keep Dry |
|  | Do Not Use if Package is Damage |  | Consult instructions for use |
|  | Do not reuse |  | Keep away from sunlight |
|  | Use By |  | Lot Number |
|  | Reference Number |  | Manufacturer |
|  | Date of Manufacture |  | Contains sufficient for <n> tests |