Sonda Evolume HER2/CEN17 FISH, kit

Evolume HER2/CEN17 FISH Probe Kit – Evolume

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Nome: Sonda Evolume HER2/CEN17 FISH, kit

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Evolume HER2/CEN17 FISH Probe Kit – Evolume

For professional use only

Product Name: Evolume HER2/CEN17 FISH Probe, kit

Packaging Specifications


Reference	Description	Kit Size
EP-13-11991	Probe HER2/CEN17 FISH, kit	5 Tests/Kit
EP-13-11992	Probe HER2/CEN17 FISH, kit	10 Tests/Kit
EP-13-11993	Probe HER2/CEN17 FISH, kit	20 Tests/Kit

INTENDED USE

The HER2 FISH Probe Kit is designed to detect amplification of the HER-2/neu gene via fluorescence in situ hybridization (FISH) in formalin-fixed paraffin-embedded (FFPE) human breast cancer tissue specimens, to provide the physician with supporting information for diagnosis.

Test Principle

Probe Design Description

The HER2/CEN17 probe is located on chromosome 17. A 400Kb probe covering the HER2 gene region is labeled with an orange dye, and part of chromosome 17 (CEN17) is labeled with a green dye. The HER2 gene marker region is located on 17q12. The result is used to guide the drug therapy and/or prognosis evaluation of primary invasive carcinoma of breast cancer.

Technique Description

Fluorescence in situ hybridization (FISH) is a molecular assay technique that uses fluorescently labeled DNA or RNA fragments or probes to hybridize via base complementary mechanism to the targeted regions of chromosomes in specimens. The probes emit fluorescence signals at a specific wavelength after excitation by a light source, and the hybridization to the chromosome can be directly observed through a microscope equipped with a proper filter set. This technique can accurately detect chromosomal abnormalities, therefore, the status of genes located at the affected regions of chromosomes can be evaluated.

Main Components

Specification	Kit content	Main components	Volume	Quantity
5 tests/kit	HER2/CEN17probe	HER2/CEN17 probe, COTIDNA, formamide,SSC,and	50 µl	1 vial
10 tests/kit	(HER2 orange probe,	dextran sulfate, etc.	100 µl	
20 tests/kit	CEN17 green probe)		200 µl	
	DAPI Mounting Medium	Antifade reagent,DAPI, and glycerin, etc.	50 µl	1 vial
			100 µl	
			200 µl	

Notes: The components in kits of different lots should not be interchanged.

Material need but not provided

- Rubber cement
- Deparaffinization solution
- Absolute ethanol
- Purified water
- 1 M HCl
- 1 N or 2 N Hidróxido de Sódio (NaOH)
- 20×SSCSolution (3 M sodium chloride, 0.3 M sodium citrate, pH 5.3)
- NP-40
- Sodiumthiocyanate
- Pepsin(250 U/mg)

Applicable instruments

- Fluorescence microscope.
- Microliter pipettor and sterile tips
- Microcentrifuge
- Water bath
- Coplin jars
- Hybridization instrument

Notes: Fluorescence microscope. The configuration of the required fluorescence microscope includes: 10× eye piece, and 10×, 40× and 100× objective lenses. It is recommended that before using the probe, the user should ask the filter set supplier for the details of the filter set to be used, so as to choose a filter set that is compatible with the labeled fluorescent dyes.

Orange fluorescence: excitation maximum: 552 nm, emission maximum: 576 nm.

Green fluorescence: excitation maximum: 496 nm, emission maximum: 520 nm.

DAPI: excitation maximum: 340 nm, emission maximum: 488 nm.

Storage Conditions and Shelf Life

Storage conditions: below-15°C, sealed, and stored in the dark.

Shipping conditions: The kits must be transported under controlled temperature conditions between 2 °C and 8 °C. The shipping time must not exceed 10 days. During transport, the temperature must not exceed room temperature.

Shelf life: 12 months.

Note: see the production date and shelf life on the outer package

Sample Requirements

1. Sample type: formalin-fixed,paraffin-embedded tissue sections, it is recommended that paraffin samples less than two years be selected to ensure the accuracy of the test.
2. Sample collection method: the fresh tissues should be fixed in 10% neutral formalin buffer for 24-48 hours within 0.5-1 hour after being removed. The tissue size should not exceed 0.5 cm³ and the recommended section thickness is 3 µm- 5 µm.
3. Sample storage: the prepared paraffin sections are recommended to be used immediately, or stored at -20±5°C.

Test Methods

A - Preparation of Working Reagen

Ethanol Wash Solutions

Prepare v/v dilutions of 70% and 85% using 100% ethanol and purified water. Dilutions may be used for 1 week unless evaporation occurs or the solution becomes diluted due to excessive use. Store at room temperature in tightly capped containers when not in use.

Pretreatment Solution

- To prepare, add together:

80 g	Sodium thiocyanate
800 mL	Purified water
1000 mL	FinalVolume

Mix thoroughly. Adjust volume to 1 liter with purified water. Filter through 0.45 µm pore filtration unit. Discard used solution at the end of each day. If the reagent is turbid or contaminated, it should be discarded immediately.

1× EDTA Solution (10 mM EDTA, 10 mM Tris-HCL)

- To prepare, add together:

20 mL	0.5 M EDTA pH8.0
5 mL	2 M Tris-HCL pH 8.0
975 mL	Purified water
1000 mL	Final Volume

Mix thoroughly. Measure pH at room temperature using a pH meter. Adjust pH to 8.1±0.1 with 1 N or 2 N NaOH. Adjust volume to 1 liter with purified water. If the reagent is turbid or contaminated, it should be discarded immediately.

2× SSC /0.3%NP-40 Washing Solution

- To prepare, add together:

100 mL	20× SSC pH 5.3
847 mL	Purified water
<u>3 mL</u>	NP-40
1000 mL	Final volume

Mix thoroughly. Measure pH at room temperature using a pH meter. Adjust pH to 7.0 ± 0.2 with 1 N or 2 N NaOH. Adjust volume to 1 liter with purified water. If the reagent is turbid or contaminated, it should be discarded immediately.

Enzyme buffer (pH 2.0 HCl)

- To prepare, add together:
- | | |
|---------------|----------------|
| 10 mL | 1 M HCl |
| <u>900 mL</u> | Purified water |
| 1000 mL | Final Volume |

Mix thoroughly. Measure pH at room temperature using a pH meter. Adjust pH to 2.0 ± 0.2 with 1 M HCl. Adjust volume to 1 liter with purified water. If the reagent is turbid or contaminated, it should be discarded immediately.

B - Preparation of instruments

Note: 2 mg/ml Pepsin solution (Take 100 mg pepsin, add 50ml enzyme buffer and mixwell) should be prepared fresh daily.

1. Turn on the thermostat water bath and set the temperature to 80°C . Pour the Pretreatment Solution into the Coplin jar, put the jar in the water bath, and preheat it to 80°C .
2. Turn on the thermostat water bath and set the temperature to 37°C . Pour 50 ml 1 mg/ml Pepsin solution into the Coplin jar, put the jar in the water bath, and preheat it to 37°C .
3. Before post-hybridization wash, set the temperature of the thermostat water bath to 68°C , pour the 2× SSC/0.3% NP-40 Washing Solution into the Coplin jar, and place the jar in the water bath for at least 30 minutes. Ensure that the temperature of the washing solution reaches $68 \pm 1^{\circ}\text{C}$ before use. Pour the 2× SSC/0.3% NP-40 Washing Solution into another Coplin jar and keep room temperature.

Testing process

Sample pretreatment

1. Preheat the slide at 65°C for 5 minutes in an oven or using the hybridization instrument or hybridizer, recommended ARA (EP-31-20147/EasyPath) or equivalent.
2. Immediately insert the preheated slide into a slide holder immersed in the deparaffinization solution (isoparaffin type) or xylene and incubate at room temperature for 15 minutes.
3. After incubation, transfer the slide holder into freshly prepared deparaffinization solution and incubate for additional 15 minutes.
4. After deparaffinization, remove the slide holder and place it on paper towel for 2 minutes to drain extra liquid.
5. Dehydrate the slide in absolute ethanol twice for 5 minutes each.
6. Rehydrate the slide in 85% ethanol, 70% ethanol and purified water in sequence for 1 minute each

7. Take out the slide and use the water absorbing paper to remove the residual solution from the slide.
8. Place the slide in the pretreatment solution preheated to 90°C and incubate for 20 minutes.
9. Take out the slide and immerse it in purified water for 1 minute. Then, take it out and remove the residual liquid.
10. Immerse the slide in 1 mg/ml Pepsin solution preheated to 37°C and incubate for 25 minutes.
11. After the enzymatic digestion, immerse the slide in purified water for 1 minute.
12. Dehydrate the slide in 70% ethanol, 85% ethanol and 100% ethanol in sequence for 1 minute each.
13. Dry the slide at room temperature.

Sample denaturation and hybridization

1. Preheat the probes to the room temperature, swirl and mix them well, centrifuge the mixture for 2s-3s in a microcentrifuge, and drop 10µl probes to the hybridization area.
2. Place the cover slide (18mmx18mm), avoiding air bubbles, and seal the slide with rubber cement.
3. Set the hybridization instrument program: denaturation: 85°C, 5 min; hybridization: 42°C, 16-18 hours.

Note: The fast hybridization program: denaturation: 85°C, 5 min; hybridization: 42°C, 2 hours.

Post-hybridization wash

1. Carefully remove the Rubber Cement.
2. Place the slide in the room temperature 2×SSC/0.3%NP-40 Washing Solution in the Coplin jar to wash off the cover slide (2min-5min).
3. Take out the slide and rinse it for 2 min with the 2×SSC/0.3%NP-40 Washing Solution preheated to 68±1°C.
4. Take out the slide and rinse it for 1 min with the 2×SSC/0.3%NP-40 Washing Solution preheated at the room temperature.
5. Take out the slide, absorb the residual solution, put the slide in 70% ethanol at room temperature, and rinse it for 1 min.
6. Air-dry the slide naturally in the dark for use later.

Counterstaining




1. Drop 10µl counterstaining solution to the hybridization area (avoiding bubbles), place the cover slide, carefully remove the excess counterstaining solution with water absorbing paper, and counter stain the slide at -20±5°C in the dark for more than 20 min.
2. Select the appropriate filter set to observe the results under the fluorescence microscope.

Slide viewing and data collection

1. Select the appropriate filter set to observe the results under the fluorescence microscope.
2. Find the area under a 10× objective lens and count the number under a 100× objective lens.
3. Adjust the focal length to find the suitable area, requiring complete nucleus boundary, uniform staining by the counterstaining solution, no overlap of nuclei, and clear signal points.

- In the selected region, find all signal points at different levels of the nuclei, and count the individual orange and green signals; the nuclei with out signals or with weak signals are not counted.

Interpretation of common signal type

	● =orange signal, HER2	● =green signal, CEN17
1		The nuclei are overlapping, so the entire areas of both nuclei cannot be seen clearly; however, the signals are not in the overlapped area, and there are two orange signals and two green signals in each nucleus
2		Two orange signals and two green signals are counted, and one of the orange signals is dispersed
3		Don't count. All nuclei overlap, so all areas cannot be seen clearly, and some signals are in the overlapped area
4		Two orange signals and two green signals are counted, and one of the orange signals is split
5		One orange signal and two green signals are counted, and one of the green signals is split
6		Two orange signals and one green signal are counted
7		Three orange signals and one green signal are counted
8		Four orange signals and two green signals are counted

Threshold determination for FISH negative specimen

Twenty specimens without the known genetic abnormalities to the FISH probe set are required to be randomly selected to prepare control slides. After hybridization, count 200 cells from each control slide and record any types of signal patterns that are considered to be FISH positive. Calculate the percentage of positive signals observed and determine the standard deviation. The threshold for FISH-negative specimens is set as the average value of percentage plus 3 times the standard deviation. The control or FISH-negative specimens are usually chosen to be the same type of tissue or cell as that of the intended target of detection by FISH.

Importance of threshold determination

The threshold must be set when probes are used for the first time since this value will serve as the reference to distinguish FISH-positive and negative specimens. If the experimental procedures are altered, such as the method of specimen pre-treatment or a change of equipment, the threshold may vary; therefore must be reset under the newly established experimental conditions.

Interpretation of Test Results

Common problems and solutions in the testing process:

Problem	Possible cause	Recommended solution
No signal or weak signal	Insufficient denaturation of samples and probes	Ensure that the temperature of the slide is $85\pm 1^{\circ}\text{C}$ during the denaturation; Extend the slide denaturation time by 2-4 min.
	No probe added	Fully thaw the probes, and ensure that the probe reagent is drawn by the pipette.
	Insufficient quantity of probes	Ensure that the probes reach room temperature before use, and ensure that the probe reagent is drawn by the pipette.
	Insufficient drying of slide	Before adding the probes onto the slide, ensure that the ethanol solution on the slide has completely volatilized.
	Too fast drying of probes	Cover the target area by the cover slide immediately after the probes are added; for the elution, remove the cover slide on only one slide at a time, and immerse the slide in the washing solution immediately before the cover slide on the next slide is removed.
	Bolhas sob a lâmina durante a hibridização	Cover the surface of the probes with the cover slide, and gently squeeze it to eliminate the air bubbles.
	Inappropriate hybridization conditions	Ensure the compliance with the specified hybridization time and temperature; do not leave any gap when sealing the slide with fluoro gel; Adjust the hybridization time and humidity as appropriate.
	Incorrect washing solution or elution conditions	Ensure that the washing solution is prepared in accordance with the IFU; ensure that the temperature of the washing solution reaches the temperature specified in the elution step; remove the cover slide before immersing the slide in the washing solution.
	Improper storage of probes or sample slides	Store the probes below -15°C in the dark; dry the unhybridized slides for long-term storage at $-20\pm 5^{\circ}\text{C}$ or short-term storage (generally no more than two weeks) at room temperature; Store the hybridized slides at $-20\pm 5^{\circ}\text{C}$ in the dark.
	Inappropriate selection of filter set for observation	Use the correct filter set to observe the probe fluorescence.
Excessively intense background of slides	Inappropriate microscope structure and objective lens for observing FISH samples, or damage of the filter set	Please contact the microscope manufacturer.
	Insufficient washing of slides before sample preparation	Immerse the slide in absolute ethanol and wipe it dry with water absorbing paper before dropping the reagent.
	Inadequate elution after the hybridization	Ensure that the washing solution is prepared correct in accordance with the IFU; operate in accordance with the IFU; ensure that the temperature of the washing solution is correct; remove the cover slide and repeat the elution step.
	Long-term use or improper storage of washing solution	Ensure that the washing solution is stored at $2-25^{\circ}\text{C}$, and discard the reagent if it is turbid or contaminated.
Too weak counterstaining	Too high or too low hybridization humidity	Adjust the hybridization humidity to the optimum.
	Too weak counterstaining	Remove the cover slide and immerse the slide in the washing solution for 5 min at room temperature. Place the slide in 70%, 85% and 100% ethanol solutions for 1 min each for gradient dehydration, and then perform counterstaining.
	Expiration or excessive light exposure of staining solution	Ensure that the staining solution is stored below -15°C in the dark; ensure that the staining solution is not expired.


Limitations of Test Methods

1. This kit is based on the fluorescence in situ hybridization technique for the HER2 gene amplification.
2. This reagent is applicable to paraffin tissue sections, and adjust the pretreatment process according to the thickness of the tissue sections or the type of samples.

Precautions

1. This kit is only for research.
2. During the operation of this kit, it is necessary to wear latex gloves to avoid the contact of the reagent with the skin. In the case of accidental contact, rinse immediately with plenty of water.
3. Samples should not be exposed to acid and alkali or be overheated, which may damage DNA and lead to failure of FISH test.
4. All components in the kit should be used within the shelf life.
5. The staining solution contains DAPI, which is a mutagen. Inhalation, ingestion or contact with the skin should be avoided.
6. The probe contains a teratogen called formamide. Contact of it with the skin and mucous membranes should be avoided.
7. In order to obtain ideal results, it is necessary to ensure that the reagents are properly prepared and stored in accordance with the IFU.

Explanation of Signs

	Upper limit of temperature		Keep Dry
	Do Not Use if Package is Damaged		Consult instructions for use
	Do not reuse		Keep away from sunlight
	Use By		Lot Number
	Reference Number		Manufacturer
	Date of Manufacture		Contains sufficient for <n> tests