#easytips

Parafina Histológica HistoPar

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1. Introdução

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User manual and best practices

1. Introduction

Paraffin is a petroleum-derived wax composed primarily of saturated hydrocarbons. In anatomical pathology laboratories, it is used in the tissue embedding process, providing support to enable thin histological sections to be made using microtomes.

However, not all paraffin is the same: currently, in addition to traditional (or common) paraffin, there are modified formulations with polymer-added paraffin that offer important benefits in terms of performance and cutting quality.

In this presentation, we will understand the differences between these two types of paraffin, their indications, and the impacts each can have on laboratory work.

2. What is paraffin used in anatomical pathology?

Paraffin is the most widely used embedding medium for histological tissue analysis. In the routine of anatomical pathology laboratories, it is indispensable for the embedding and cutting of samples, enabling the production of slides with quality and diagnostic precision.

Tissue embedding begins after proper fixation and dehydration.

Paraffin's insolubility in water means that specimens need to be completely dehydrated through increasing concentrations of alcohol (70%, 90%, and 100%). The process requires clearing agents such as Xylene. This vital step makes tissues translucent and ready for complete paraffin infiltration.

During processing, paraffin in its liquid state infiltrates the interior of the tissue and, upon rapidly solidifying during cooling, forms a firm matrix. This matrix surrounds and stabilizes the sample, preserving the architecture and preventing distortions of the tissue structure during the microtomy stage.

To be used in histology, paraffin must meet the following requirements:

- ▶ Appropriate melting temperature (generally between 54°C and 60°C)
- ▶ Low impurity content
- ▶ Capability for homogeneous tissue infiltration

Heating paraffin to an elevated temperature alters the paraffin's properties.

3. Common Paraffin (Without Polymer)

Main Characteristics

Traditional paraffin is basically composed of pure paraffinic hydrocarbons. Among its characteristics are:

- ▶ Good plasticity
- ▶ Easy melting and rapid solidification
- ▶ Good compatibility with histological stains

However, it may present limitations such as:

- ▶ Formation of cracks or breaks during cutting, especially in more fragile tissues
- ▶ Greater difficulty in obtaining ultrathin, serial, or cracked sections
- ▶ Lower elasticity, which can generate folds
- ▶ Batch inconsistency
- ▶ Greater wear on the microtome blade edge, increasing costs
- ▶ Higher melting point between 59°C and 62.8°C, requiring increased temperature in the processor and embedding center for melting to occur

Common Paraffin (Without Polymer)

4. Polymer Paraffin: What Changes?

Polymer paraffin is the result of adding synthetic resins or special polymers to pure paraffin. This modification aims to improve important physical properties, such as:

- ▶ Elasticity
- ▶ Cohesion
- ▶ Tissue support capacity during cutting

Observed Benefits:

- ▶ Greater resistance to breaks and cracks
- ▶ Ease in obtaining continuous serial sections
- ▶ Reduction of artifacts (folds, cracks, disaggregations)
- ▶ Improved cutting quality, especially in fragile, small tissues or those with high fat content
- ▶ Batch consistency, as it is formulated exclusively for histology
- ▶ Reduced wear on the microtome blade edge, reducing costs
- ▶ Low melting point of 56°C to 60°C
- ▶ Change in block coloration if temperature exceeds 65°C to yellowish coloration, ensuring that the sample is processed and embedded at ideal temperature, without damage to the sample

Paraffin with added polymer

Figure A: At temperatures up to 65°C, the paraffin coloration does not change, ensuring sample preservation.

Figure B: Yellowish due to temperature above 65°C. The higher the temperature, the greater the intensity of the yellowish color, indicating quality deviation in sample processing and embedding.

- The higher the temperature, the more pronounced the yellowish tint will be, indicating a deviation in quality during processing and inclusion of the sample.

4.1 Practical tips for best use

- ▶ **Tissue processor bowl temperature:** up to 65°C, preferably 62°C
- ▶ **Embedding Center Temperature:** up to 65°C, preferably 62°C
- ▶ Use the larger wheel, making faster rotations to facilitate formation and maintain the integrity of the histological ribbon
- ▶ **Histological Bath:** up to 40°C, as temperatures above this can cause "fraying" or melting, impairing microscopy

5. Technical Guidance – Equipment Calibration and Checking

Practical Tips for Best Use

To ensure the correct use of **Histopar** polymer paraffin from **EasyPath**, ensuring the preservation of its physicochemical properties and obtaining high-quality histological results, it is necessary to calibrate your equipment correctly.

Histopar should be used at a **maximum temperature of 65°C**, and it is essential that the equipment involved in the histological process be properly adjusted and calibrated to operate within this limit.

Equipment and Recommended Temperatures

Equipment	Maximum Operating Temperature
Paraffin processor	Up to 65 °C
Embedding center	Up to 65 °C
Histological bath	Up to 40 °C

Importance of Temperature Checking

Many laboratories currently use paraffins with higher melting points, generally without addition of plastic polymers and not specifically developed for histological use. These paraffins require higher temperatures (above 65°C) for complete melting, as their melting point is above 62°C. When the same parameters are maintained for Histopar, overheating occurs, compromising its polymer structure and resulting in loss of **elasticity** — an essential property for good handling and cutting of blocks.


Figure C 40°C - Histological Bath

Figure D: 65°C - Liquid Paraffin

6. Problems Identified in Laboratories

6.1 Use of Inadequate Temperatures

For paraffins without polymers, equipment is often adjusted above 70°C or even 80°C to ensure complete melting. However, these temperatures are **incompatible** with Histopar's composition and can **damage tissue structure and alter paraffin properties**.

6.2 Incorrect Equipment Calibration

It is common for the **digital or analog display to indicate 65°C** while the actual paraffin temperature is significantly higher. The absence of periodic calibration leads to inaccurate readings, causing inadvertent use of temperatures above the recommended limit.

6.3 Excessive Heating to Accelerate the Process

In some cases, users purposely raise the temperature to accelerate paraffin melting. However, temperatures above 65°C **cause irreversible alteration of the polymer structure**. Even if the temperature is subsequently reduced, **Histopar's elasticity will already be compromised**, impairing its application.

7. Operational Recommendations

- Verify and adjust the actual temperature of equipment before use, using a calibrated external thermometer
- Regularly calibrate the tissue processor, embedding center, and histological bath according to manufacturer's guidance
- Maintain daily temperature control through recording in check sheets
- Avoid intentional overheating with the objective of accelerating the paraffin melting process
- Report temperature deviations above 65°C to the responsible technical team before starting to use Histopar

8. Conclusion

The proper use of **Histopar paraffin** requires **special attention to temperature control of histological equipment**.

Observance of maximum recommended temperatures — **65°C for the processor and embedding center, and 40°C for the histological bath** — is **crucial** to preserve polymer properties and ensure the integrity of processed tissues.

Failure to comply with these guidelines may lead to loss of sample quality and false perception of product failure when, in reality, the problem stems from **incorrect use of temperature parameters**.

9. Frequently Asked Questions

Does polymer paraffin alter histological staining?

No if it is used correctly and renewed regularly.

Can I mix common paraffin and polymer paraffin?

It is not recommended, as it may compromise block homogeneity.

Does the microtome need adjustment?

In some cases, yes. It is recommended to use the larger wheel, making faster rotations to facilitate formation and maintain the integrity of the histological ribbon.

Does any product need to be added to the paraffin?

No, polymer paraffin is ready to use, in addition to being filtered and purified, ensuring better quality for the examination.

Why does polymer paraffin change color at temperatures above 65°C?

This is an additional and easily visible protection to ensure sample integrity, as above 65°C sample deterioration begins, which is observed mainly in more complex examinations such as immunohistochemistry or molecular biology.

10. Bibliographic References

Lim JI, Lim KJ, Choi JY, Lee YK. Modified paraffin wax for improvement of histological analysis efficiency. *Microsc Res Tech*. 2010 Aug;73(8):761-5. doi: 10.1002/jemt.20818. PMID: 20063315.

Suvarna SK, Layton C, Bancroft D. Bancroft's theory and practice of histological techniques. 7th ed. 2019. doi.org/10.1016/C2015-0-00143-5

This material is the property of Grupo Erviegas. Its total or partial reproduction, by any means, is not authorized without prior authorization from Grupo Erviegas.

Produced by the Research and Development department - Grupo Erviegas

Erviegas Química Fina e Plásticos LTDA — CNPJ 46.271.011/0001-07
Al. Plutão, 593, American Park Empresarial — Indaiatuba/SP, Brazil — ZIP: 013347-656
Phone: (11) 5034-2227
Website: www.grupoerviegas.com.br
Help and Support: ajuda.grupoerviegas.com.br